Oxyguard Reusable

Instructions for Use

OXYGUARD®

Endoscope Biteblock with Oxygen Access

INTENDED USE: The Oxyguard is a biteblock intended to protect the patient and instrument during endoscopy. It includes an oxygen port to allow supplemental oxygen to be delivered simultaneously to the nose and mouth.

The Oxyguard is intended for use by medical staff competent in the management and maintenance of a clear airway in an unconscious or sedated person. It is intended for adult patients.

INSTRUCTIONS FOR USE

1. Remove the biteblock from the packaging.
2. If required, attach a head strap to the notches on either side of the biteblock.
Note: Use of a head strap is optional. It may help keep the biteblock in position during endoscopy.
3. If supplemental oxygen is required, attach the oxygen line to the port on the side of the biteblock (see Figure 1).
4. Position the Oxyguard in the patient's mouth (see Figure 2), and fit and adjust the head strap (if one is used).
5. Adjust the oxygen flow accordingly.
Note: Oxygen will flow to the patient's nose and mouth.
6. If a head strap is used, detach this from the biteblock first before removing the biteblock from the patient's mouth.

WARNING: The Oxyguard should only be re-used between patients after it has been appropriately high-level disinfected or sterilized. This is to avoid the risk of cross-infection.

AUTOCLAVING

The Oxyguard can withstand up to 12 autoclave cycles before the device may start to weaken.

Before each use, visually examine to check the ongoing integrity of the device.



Prior to autoclaving, the device must be thoroughly manually cleaned (eg, with a brush and detergent), paying particular attention to channels and crevices. All traces of dirt and/or secretions must be removed prior to autoclaving.

Note: When ready for disposal, ensure the product is appropriately disposed of as a contaminated biohazard.

WARNING: Consult your infection control department or infection control expert prior to autoclaving and re-use of this device.

WARNING: The black head strap contains natural rubber latex. Latex has been known to cause allergic reactions.

CAUTION: US federal law restricts sale of this device to, or on the order of, a physician.

Trawax Pty Limited
28 Queenscliff Road
Queenscliff
NSW 2096, Australia
www.trawax.com.au

Distributed in Australia by:
Vertrieben in Australien durch/Distribué en Australie
par :

Endomed Pty Ltd
Ph: +61 7 3881 1883
Fx: +61 7 3881 0211

Advena Ltd.
Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

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©Trawax Pty Ltd 2007. Rev 09 11/19
Oxyguard is a registered trademark of Trawax Pty Ltd.